

Phone: 866.592.0488 • Fax: 209.834.8832
Tracy, CA 95304
1851 East Paradise Road, Suite A



Deluxe Tendonitis Strap

Instructions For Use

Please read these instructions completely and carefully before applying this brace. Correct application of the brace is important for proper function.




Deluxe Tendonitis Strap

Instructions For Use

Please read these instructions completely and carefully before applying this brace. Correct application of the brace is important for proper function.

Indications:

Medial epicondylitis (golfer's elbow) and lateral epicondylitis (tennis elbow).



1. Wrap support around forearm at least two finger widths below crease of elbow.
2. Loop end through D-Ring and secure with pad on the affected side of the arm.
3. Support should fit snug, but not tight enough to impair circulation.

Use and Care: Hand wash in cold water using mild soap, rinse thoroughly and air dry. If not rinsed thoroughly residual soap may cause skin irritation and deterioration of the material.

Warning: This product is to be used under the supervision of a medical professional. This device is not intended for unsupervised public use. If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

Warning: Neoprene supports should not be worn by individuals with known susceptibility to dermatitis or allergies to neoprene.

Warranty: Top Shelf Orthopedics guarantees this product to be free of defects in material and construction for a period of 6 months from the date of purchase.

Caution: Federal Law (U.S.A.) restricts this device for sale by or on the order of a licensed healthcare professional.