

topshelforthopedics.com  
Phone: 866.592.0488 • Fax: 209.834.8832  
Tracy, CA 95304  
1851 East Paradise Road, Suite A



## Horizon Walker Instructions for Use

Please read these instructions completely and carefully before applying this product. Correct application is important for proper function of this brace.



# Horizon Walker

## Instructions for Use

Please read these instructions completely and carefully before applying this product. Correct application is important for proper function of this brace.

**Indications:** Soft tissue injuries of the lower leg, acute ankle sprains, stress fractures of the lower leg, stable fractures of the foot and ankle, stable delayed union or non-union fractures of distal tibia and fibula, achilles tendon repairs.



1. Open the boot liner without removing it from the boot.
2. Place the boot flat on the floor and put the foot inside the boot. Make sure the heel is touching the back of the boot.
3. Secure the lining to the leg by overlapping the liner on the top of the foot and at the front of the lower leg.
4. Secure the Velcro straps beginning with the straps closest to the toes.

### **Use and Care:**

Hand wash the liner with mild soap and water and lay out to dry. Do NOT machine dry.

### **Notes:**

For single patient use only.

### **Warning:**

This product is to be used under the supervision of a medical professional. If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

### **Warranty:**

Top Shelf Orthopedics guarantees this product to be free of defects in material and construction for a period of six months from the date of purchase.

### **Caution:**

Federal Law (U.S.A.) restricts this device for sale by or on the order of a licensed healthcare professional.