

Forefoot Offloading Shoe

The forefoot offloading shoe was designed to effectively reduce metatarsal head pressure and great toe pressure after forefoot injuries, surgery, or trauma.



Indications

- Post-surgical healing
- Forefoot trauma
- Metatarsal head & phalangeal decompression
- Forefoot offloading

Features & Benefits

- Wedge sole reduces forefoot pressure by shifting body weight to the midfoot and heel
- Square toe design to protect the toes
- Ankle strap holds foot firmly in the shoe to reduce motion of the foot within the shoe
- Adjustable forefoot closure can comfortably accommodate bandages

Part #	Size	Men's Shoe Size	Women's Shoe Size
501762	X - Small	N/A	4 - 7
501763	Small	6 - 8	7.5 - 10
501764	Medium	8.5 - 10	10.5 - 13
501765	Large	10.5 - 12	N/A
501766	X - Large	12.5 - 14	N/A



Ph: 866.592.0488

Fax: 209.834.8832

1851 E. Paradise Rd., Ste A

Tracy, CA 95304

To view our full product line, please visit us at
www.topshelforthopedics.com