

Intelli-Flo Hip Pad

Joint Specific

Patented in-pad flow restrictor tailors cold temperature specifically for the hip

Optimal Flow


Flow channels direct the cold for optimal water circulation throughout the pad

Targeted Cooling Relief

Customized pad provides additional cold where it's needed most - at the incision


Comfortable

Ergonomic pad design contours the hip for optimal comfort and compliance



Intelli-Flo Technology

The Intelli-Flo hip pad provides custom-tailored temperature control specific to the hip. The pad is designed to keep a safe operating temperature range within 45°F - 55°F.



Intelli-Flo Hip Pad Indications for Use:

- Post Operative
 - Total Hip Procedures
 - Partial Hip Procedures
 - Hip Arthroscopy Procedures
- Post-Trauma
- Chronic Pain
- Physical Therapy



Distributed by:
 Pacific Medical, Inc.
 1700 N. Chrisman Road
 Tracy, CA 95304
 Phone: 866.592.0488 Fax: 209.834.8832

Improves Patient Recovery

Polar Care Kodiak cold therapy provides pain management throughout post-operative care and rehabilitation, and helps patients increase range of motion during the recovery process. Unlike messy ice, the Polar Care Kodiak is a closed system that transfers water from the unit to the joint specific Intelli-Flo pads minimizing the risk of infection.

POLAR CARE KODIAK

Part No.	Description
502700	Polar Care Kodiak

POLAR CARE KODIAK INTELLI-FLO HIP PAD

Part No.	Description
502721	Intelli Flo Pad, Hip
502722	Intelli Flo, Hip, Sterile Dressing



Top Shelf Orthopedics
 Tracy, CA 95304
 Phone: 866.592.0488
 Fax: 209.834.8832

