



 TELE-RANGE

 TOP SHELF
ORTHOPEDICS™

Post-Op Bracing

The Tele-Range post-operative brace family provides an unparalleled design, which accommodates a variety of patients, and protocols.



Indications:

- ACL, PCL, MCL, LCL surgeries and injuries
- Meniscal Repairs
- Osteo Chondral Transfers and/or repairs
- High Tibial Osteotomy
- Patella Tendon repairs
- Stable Tibial Plateau Fracture

Features:

- Visually interactive hinge-allows easy hinge adjustment.
- Range of motion hinge-control from $\sim 10^\circ$ to 120° in 10° increments.
- Quick lock adjustment-allows the brace to be locked with one button from 0° to 40° in 10° increments.
- Telescoping hinge bars-easy to adjust and requires no tools.
- Popliteal strap adjustment-provides accommodation for post-op dressing and swelling.
- Available in full foam version.

Part Number	Descriptions
505501	Tele-Range, Post-Op Brace
505502	Tele-Range XL, Post-Op Brace
505503	Tele-Range, Full Foam, Post-Op Brace
505504	Tele-Range, Full Foam, XL, Post-Op Brace

Post-Operative Bracing Accessories

505505	Tele-Range Replacement Pad Set
505506	Tele-Range, Full Foam Replacement Pad Set
505507	Tele-Range, Full Foam XL Replacement Pad Set
505508	Tele-Range, Replacement Condyle Pad Set

