

TR Post-Op Brace Family

Featuring the Tele-Range hinge, the upper and lower hinge bars can be easily shortened or contoured for optimum patient fit and comfort.

Indications:

- ACL, PCL, MCL, LCL injuries
- Tibial plateau fractures
- Osteochondral repairs
- Meniscal repairs
- Patella tendon repairs
- Condylar fractures
- Sprains/strains of the knee
- Total knee replacement
- HTO (high tibial osteotomy)

TR Post-Op Brace

- Reduced foam is ideal for warmer climates
- Range of motion control from $\sim 10^\circ$ to 120° in 10° increments
- Quick lock adjustment allows brace to lock with one button from 0° - 40° in 10° increments
- Scored anodized aluminum hinge bars for easy breakdown to rehab style brace
- Latex free materials



Part Number	Size	Description
505511	Fits thigh up to 27"	TR, Post-Op Brace
505512	Fits thighs 27" - 35"	TR, Post-Op Brace, XL

TR Post-Op Full Foam

- Full foam for compression and soft tissue containment
- Range of motion control from $\sim 10^\circ$ to 120° in 10° increments
- Quick lock adjustment allows brace to lock with one button from 0° to 40° in 10° increments
- Scored anodized aluminum hinge bars for easy breakdown to rehab style brace
- Available in XL version
- Latex free materials

Part Number	Size	Description
505513	Fits thigh up to 27"	TR, Post-Op Brace, Full Foam
505514	Fits thighs 27" - 35"	TR, Post-Op Brace, Full Foam, XL



Post-Operative Bracing Accessories

505515	TR, Post-Op Replacement Pad Set
505516	TR, Post-Op Replacement Pad Set, XL
505517	TR, Post-Op Full Foam Replacement Pad Set
505518	TR, Post-Op Full Foam Replacement Pad Set, XL

Rapid Fit Post-Op

A sized easy to apply post-op brace featuring the Tele-Range hinge.

- Range of motion control $\bar{10}^{\circ}$ to 120° in 10° increments
- Quick lock adjustment allows brace to lock with one button from 0° to 40° in 10° increments
- Full $\frac{1}{2}$ " foam provides comfort, compression and soft tissue containment
- Anodized hinge bars
- Pre-sized Small through X-Large for easier application
- Four straps with Quick fit buckles for easy on/off
- Latex-free

Sizing

Measure circumference of thigh 6" (15cm) above mid patella.

Part Number	Size	Description
505523	14" - 18"	TR, Rapid Fit Post-Op, SM
505524	18" - 22"	TR, Rapid Fit Post-Op, MD
505525	22" - 27"	TR, Rapid Fit Post-Op, LG
505526	27" - 31"	TR, Rapid Fit Post-Op, XL

